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The second generation of the range of reusable GENIUS II Shaver Blades for minimally invasive surgery - for example on knee, hip and shoulder bones - consists of completely redesigned products which promote the success of the operation more effectively than ever:

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Improved colour coding for easier recognition and allocation of the components - even in the darkened operating theatre.

- | | | |
|-----------|-----------|-----------|
| ■ GE 0135 | ■ GE 0142 | ■ GE 0155 |
| ■ GE 0242 | | |
| ■ GE 0335 | ■ GE 0342 | ■ GE 0355 |
| ■ GE 0435 | ■ GE 0442 | ■ GE 0455 |
| ■ GE 0542 | | |
| ■ GE 0642 | ■ GE 0655 | ■ GE 0665 |
| ■ GE 0755 | ■ GE 0765 | |
| ■ GE 0855 | | |

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Minimum clogging thanks to the redesigned opening at the end of the outer tube.

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Even less vibration in the motor thanks to calibrated tubes.

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Excellent cutting ability thanks to the advanced manufacturing technology and the top-quality material used, in combination with the improved blade geometry of the bevel edge cutters and the tapered/oval burs.

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- Bevel edge cutter:
More aggressive thanks to new blade configuration.
- Tapered/oval bur:
Cross-cut for fragmenting chips.

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The reverse of the outer tube is marked for improved orientation in the operative site.

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
Compatible with Linvatec, Arthrex, FMX, Dyonics EP1, Aesculap and Stryker 12K connections.

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
It's their special three-component design that makes the KOMET MEDICAL GENIUS II Shaver Blades so unmistakable. Thanks to this unique feature, the inner tube of all shavers can be cleaned with a plastic brush. Furthermore, this distinctive construction allows a visual inspection of the inner tube.

Contamination can be spotted before the sterilization procedure, in which case the instrument can be subjected to another cleaning cycle. Without this feature, neither manual cleaning nor control would be possible and the result of the preparation would therefore be uncertain.

9835 M
带有固定夹具的灭菌盒
Steri-box with clamping device
203 x 114 x 30 mm

* 套管外径
External diameter
of the outer tube

直径 mm Diameter mm	与以下接口匹配: compatible with the following connections:	Linvatec (Arthrex, FMS)	Dyonics	Aesculap	Stryker	转速 rpm Speed	刀片类型 Blade Type
3,5		GE0135.LI	GE0135.DY	GE0135.AC	GE0135.ST	6000	01 = 带齿半月板刨刀 Aggressive Meniscus Cutter
4,2		GE0142.LI	GE0142.DY	GE0142.AC	GE0142.ST	6000	02 = 双刨刀 BiCut
5,5		GE0155.LI	GE0155.DY	GE0155.AC	GE0155.ST	6000	03 = 全径刨刀 Full Radius Resector
4,2		GE0242.LI	GE0242.DY			6000	04 = 全径带齿刨刀 Aggressive Full Radius Resector
3,5		GE0335.LI	GE0335.DY	GE0335.AC	GE0335.ST	6000	05 = 斜边刨刀 Bevel Edge Cutter Resector
4,2		GE0342.LI	GE0342.DY	GE0342.AC	GE0342.ST	6000	06 = 圆头打磨器 · Round bur
5,5		GE0355.LI	GE0355.DY	GE0355.AC	GE0355.ST	8000	07 = 圆桶形打磨器 · Oval bur
3,5		GE0435.LI	GE0435.DY	GE0435.AC	GE0435.ST	8000	08 = 锥形打磨器 · Tapered bur
4,2		GE0442.LI	GE0442.DY	GE0442.AC	GE0442.ST		
5,5		GE0455.LI	GE0455.DY	GE0455.AC	GE0455.ST		
4,2			GE0542.DY		GE0542.ST		
4,2		GE0642.LI	GE0642.DY	GE0642.AC	GE0642.ST		
5,5		GE0655.LI	GE0655.DY	GE0655.AC	GE0655.ST		
6,5		GE0665.LI	GE0665.DY				
5,5		GE0755.LI	GE0755.DY	GE0755.AC	GE0755.ST		
6,5		GE0765.LI	GE0765.DY				
5,5		GE0855.LI	GE0855.DY	GE0855.AC	GE0855.ST		

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Linovatec is a trademark of ConMed Corporation. Dyonics is a trademark of Smith & Nephew Inc. Aesculap is a trademark of Aesculap AG & Co. KG. Stryker is a trademark of Stryker Corporation.
KOMET MEDICAL is not affiliated with any of the above.

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Sterilization and maintenance

*The Shaver Blades are perfectly suited for
their specific surgical use. They are designed,
constructed and produced with utmost care.*

*The Shaver Blades are suitable for reuse.
These quality instruments assure the best
working results and a long service life, provided
they are used in the proper manner.
Special care must therefore be taken during
their use and preparation.*

*Improper use may damage tissue and may lead
to premature wear, destruction of the instru-
ments and injury to the operator, patient or
other persons. The below instructions for clean-
ing and use have to be closely observed.*

***A validated method can be downloaded
from our website:***

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GENIUSadapter Validated cleaning

The newly developed GENIUSadapters provided by KOMET MEDICAL are an auxiliary device to ensure **validated cleaning**, making the cleaning result of GENIUS II Shaver Blades, in conjunction with a validated cleaning method, repeatable at any time. This means that our GENIUS II Shaver Blades - designed for reprocessing and therefore without hard-to-access areas - can be prepared in total safety. This was validated and confirmed by co. SMP GmbH.

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Instructions for use and recommendations for cleaning Shaver Blades

Indications:

Depending on their nominal diameter, the reusable Shaver Blades can be employed in various joints (e.g. knee, hip, shoulder, ankle and elbow joints and wrists). They have been developed especially for the resection and treatment of tissue (for example bones, menisci, tendons, ligaments etc.).

Contraindications:

All the general contraindications of arthroscopic joint surgery apply.

Suitable power systems:

GENIUS II Shaver Blades are compatible with the following power systems: Dynonics EP-1, Stryker 12K, Aesculap, Linvatec, Arthrex, FMS and Concept. When used properly there is no risk of damage to the respective power system.

The legal provisions governing warranties under the sales law remain unaffected by the use of compatible instruments, as long as the manufacturing tolerances are adhered to as precisely as with the original product. This fact was confirmed by the law firm Krieger Mes & Graf von der Groeben (Düsseldorf). Expert opinion available on request.

Materials:

The shavers are composed of adaptors made of blue anodised aluminium which allow anchoring in the hand piece, as well as metal tubes with stainless instrument steel blades at the inner and outer tube. The existing o-rings are made of rubber.

Cleaning and disinfecting:

- Prior to first use, new GENIUS II Shaver Blades have to undergo the same cleaning and disinfecting procedure as used instruments (see preparation cycle).
- In order to prevent residues from drying in the narrow lumen, clean and disinfect Shaver Blades immediately after use.
- The individual components must be rinsed under a water jet or in a suitable cleaning and disinfecting solution until there is no more visible contamination.
- We recommend a preparation of the Genius II Shaver Blades in the thermo disinfectant with a cleaning agent with proven material compatibility. We recommend deconex 28 Alka one (made by co. Borer Chemie).
- In order to avoid damage to the material caused by lime stains, fully demineralised water should be used for the final rinse.
- The GENIUS II Shaver Blades are to be cleaned with medically clean compressed air.
- Make sure not to mistake the blades.
- After cleaning, the components are to be examined visually to ensure that they are clean, ready for use and in perfect condition.
- Replace damaged, bent, worn or blunt Shaver Blades with new GENIUS II Shaver Blades. Repeat cleaning process if contamination persists. Do not under any circumstances exchange individual components only (for example just the inner blade!)
- After completed examination of models GEXXXX.LI, GEXXXX.DY and GEXXXX.ST, the blue adaptor and the inner tube are remounted.

Prior to sterilization:

- The Shaver Blades must be completely clean prior to sterilization.
- Lubricate the inner tube with a suitable oil-based lubricant in order to guarantee low friction.
- Place the components in a suitable container. We recommend sterilizing the Shaver in disassembled condition.

Sterilization:

- The GENIUS II Shaver Blades are now ready to be sterilised according to a validated, approved method. The recommended method is vapour sterilization in the autoclave. Not suitable for sterilization by hot air and plasma sterilization. Sterilization at low temperature should only be carried out by certified/validated preparation services.

After sterilization:

- The Shaver Blades should ideally be stored in a sealed, sterile container.
- Prior to the next surgical use, the inner tube is inserted back into the corresponding outer tube and locked by applying axial pressure (in model GEXXXX.ST locking is not possible). Make sure that only components with identical colour codes are combined.
- The GENIUS II Shaver Blades are now ready for use in your Shaver power system. For further instructions please refer to the instructions for use provided with your Shaver power system.

Precautions to be taken / Warnings:

- Never allow GENIUS II Shaver Blades to run dry, i.e. without rinsing agent, as this may lead to jamming and ensuing destruction of the inner blade.
- The length of surgical use is determined by the type and amount of tissue to be removed. Immediately exchange the GENIUS II Shaver Blade if there is obvious blunting of the blades or in case of poor performance.
- A fading of its blue colour does not indicate a deterioration of the quality of the Shaver adaptor. This merely cosmetic change is caused by unsuitable cleaning agents.
- For safety reasons it is imperative to observe the indicated speed. The permitted maximum speed should only be used for short periods of time. High speed decreases the cutting efficiency, increases the tendency to clog and leads to reduced durability of the instrument.
- Non-observance of the recommended speed for GENIUS II Shaver Blades may lead to attrition of the metal. Excessive speed also increases the risk of the instrument seizing.

For additional information regarding the proper use and preparation of the instruments please refer to the brochure "Proper Maintenance of Instruments" issued by the working group "Preparation of instruments". (AKI - www.a-k-i.org)

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Brasseler DMI (Shanghai) Ltd.