























编号系统 · 钻头

Numbering system · Rotary instruments



材质 · Material	形状 · Shape	分型 · Versions	杆部 · Shank
NC = 钨碳合金 Tungsten carbide	01 = 圆形 · Round 03 = 圆柱形 · Cylinder	例 · Examples 01 = Ø 1,0 mm 6 刻刃数 6 blades	001 = L 37,0 mm 002 = L 44,5 mm 003 = L 48,0 mm 004 = L 51,0 mm 006 = L 57,0 mm 007 = L 62,0 mm 008 = L 67,0 mm 009 = L 70,0 mm 010 = L 73,0 mm 011 = L 80,0 mm 012 = L 95,0 mm 026 = L 125,0mm 027 = L 65,0 mm 036 = L 75,0 mm 039 = L 44,5 mm 014 = L 44,5 mm 015 = L 48,0 mm 016 = L 55,0 mm 017 = L 62,0 mm 018 = L 65,0 mm 019 = L 67,0 mm 022 = L 75,0 mm 023 = L 70,0 mm 060 = L 44,0 mm 061 = L 51,0 mm 062 = L 64,0 mm 503 = L 75,0 mm 029 = L 44,5 mm
NS = 不锈钢 Stainless Steel	04 = 圆柱形 圆头 Cylinder round 04 = 橄榄形 · Olive	02 = Ø 1,4 mm 6 刻刃数 6 blades 07 = Ø 7,5 mm 8 刻刃数 8 blades 08 = Ø 9,0 mm 6 刻刃数 6 blades	Rund · Round *1
ND = 金刚砂 Diamond	06 = 细长形 · Tapered 07 = 细长形 圆头 Tapered round 08 = 细长形 尖头 Tapered pointed 09 = 蛋形 · Egg 10 = 火焰形 · Flame 11 = 火焰形 · Flame 12 = 梨形 · Pear 14 = Shannon 16 = Swanson 17 = 穿线钻 18 = Lindemann 19 = 柱形 皮肤打磨. Dermabrader cylind. 20 = 椎形 皮肤打磨 Dermabrader tapered 21 = 颅骨铣刀. Cranio 22 = 蓓蕾形 圆头 Bud round 23 = Lindemann 切割形 Lindemann Pin cutter 24 = 桶形-Barrel 25 = 橡子形-Acorn 26 = 金刚石切割盘 Diamond disc 27 = 扭纹钻 Twist Drill 27 = 螺纹钻 Thread cutter 32 = 环形切割钻 Trepan burs	88 = Ø 4,0 mm 交叉刃口 with cross-cut	Stryker *2 Xomed Dental *3 Osteon/Imperator Aesculap Hilan Stryker TPS *4
		*1 = Storz; Aesculap; Micro Aire; 3 M; Bien Air; ... *2 = De Puy/Stryker *3 = Synthes *4 = Hall Ultra Power, Midas Rex 1/2, Stryker TPS	070 = L 30,0 mm 071 = L 38,0 mm 072 = L 48,0 mm 075 = L 48,0 mm 099 = L 48,0mm* *Imperator 511 = ▲ 512 = I 513 = II 514 = III 515 = 0 516 = M 517 = AM

不锈钢 · 圆头
Stainless Steel · Round

NS 01...	直径 mm	叶片数	头部长度 mm	长度mm	长度mm	长度mm	长度mm	长度mm	长度mm
				57	70	95	125	55	44,5
				圆柄*				Stryker	Synthes
	0,6	6			NS0148.009	NS0148.012	NS0148.026		
	0,7	6			NS01100.009				
	0,8	6			NS0149.009	NS0149.012	NS0149.026		
	1,0	6			NS0101.009	NS0101.012	NS0101.026		
	1,4	8		NS0150.006	NS0150.009	NS0150.012	NS0150.026		
	1,6	6					NS0103.026	NS0103.016	
	1,8	8		NS0151.006	NS0151.009	NS0151.012	NS0151.026		
	1,8	18						NS0113.016	
	2,3	8		NS0115.006	NS0115.009	NS0115.012	NS0115.026	NS0115.016	
	2,7	10		NS0152.006	NS0152.009	NS0152.012	NS0152.026		
	3,1	6			NS0107.009		NS0107.026	NS0107.016	NS0107.029
	3,1	10			NS0153.009	NS0153.012	NS0153.026		
	3,5	12			NS0154.009	NS0154.012	NS0154.026		
	3,5	6			NS0191.009		NS0191.026		
	4,0	12		NS0120.006	NS0120.009	NS0120.012	NS0120.026	NS0120.016	
	4,5	6			NS0121.009		NS0121.026	NS0121.016	
	4,5	12			NS0123.009	NS0123.012	NS0123.026	NS0123.016	
	5,0	8			NS0109.009		NS0109.026	NS0109.016	
	5,0	12			NS0155.009	NS0155.012	NS0155.026		NS0155.029
	6,0	14			NS0156.009	NS0156.012	NS0156.026	NS0156.016	
	7,0	16			NS0125.009	NS0125.012	NS0125.026	NS0125.016	

* = Storz; Aesculap; Micro Aire; 3 M; Bien Air;...