

杆部类型 · Shank type

圆柄 · Round shank Ø 2.35 mm

总长 · Overall length

002		44.5 mm
004		51.0 mm
006		57.0 mm
009		70.0 mm
012		95.0 mm
026		125.0 mm
001		37.0 mm
003		48.0 mm
007		62.0 mm
039		44.5 mm

齿科器械杆部类型 · Dental shank type

029		44.5 mm
-----	--	---------

Imperator杆部 · Imperator shank

099		70.0 mm
-----	--	---------


CRANIO

025		44.5 mm
037		51.0 mm

XOMED 杆部 · XOMED shank

Stryker杆部 · Stryker shank

(J-notch)



总长 · Overall length

060	44.0 mm
061	51.0 mm
062	64.0 mm
503	75.0 mm



总长 · Overall length

014	44.5 mm
015	48.0 mm
016	55.0 mm
017	62.0 mm
018	65.0 mm
019	67.0 mm

高速旋转器械 · High-Speed

Stryker杆部 · Stryker shank (Old Midas Rex, Hall Ultra Power)



515	O
516	M
517	AM

